

Product Datasheet

NED-3238 (orb1910833)

Description	NED-3238 is highly potent third generation inhibitors of human arginase I and II with IC50 of 1.3 nM...
Note	For research use only
Formula	$C_{17}H_{30}BCl_2N_3O_4$
MW	422.154
SMILES	<chem>O=C([C@]1(N)CN(C[C@@H](N)CC2=CC=CC=C2)C[C@@H]1CCCB(O)O)O.[H]Cl.[H]Cl</chem>
Expiration Date	12 months from date of receipt.