

Product Datasheet

L343 (orb1911045)

Description	L343, an ionizable cationic lipidoid, facilitates the synthesis of liposomes for the systemic delivery of RNAi therapeutics.
Storage	-20°C
Note	For research use only
Formula	$C_{51}H_{95}NO_6$
Purity	98.00%
MW	818.3
CAS Number	1443520-22-6
SMILES	<chem>CN(C)CCCC(=O)OC(CCCCCCCC\C=C/CCCCCCCC(=O)OC(C)(C)CCCCCCCC\C=C/CCCCCCCC(=O)OC(C)(C)C</chem>
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)