



Product Datasheet

TCL1A antibody (orb191456)


Description	Rabbit polyclonal antibody to TCL1A
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human TCL1A(NP_001092195.1).
Preservatives	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	16kDa
Uniprot ID	P56279
Dilution Range	WB: 1:500-1:2000, IHC-P: 1:50-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of Ramos cell lysa...



Immunohistochemical analysis of formalin...



Immunohistochemical analysis of formalin...