

Product Datasheet

hnRPD (Phospho-Ser83) Antibody (orb191551)

Description	hnRPD (Phospho-Ser83) Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	Q14103
Expiration Date	12 months from date of receipt.