

Product Datasheet

IL-2R β /CD122 (Phospho-Tyr364) Antibody (orb191552)

Description	IL-2R β /CD122 (Phospho-Tyr364) Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	P14784
Expiration Date	12 months from date of receipt.