

Product Datasheet

MER/SKY (Phospho-Tyr749/681) Antibody (orb191589)

Description	MER/SKY (Phospho-Tyr749/681) Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	Q12866
Expiration Date	12 months from date of receipt.