


Product Datasheet

Trk B (Phospho-Tyr706/Tyr707) Antibody (orb191605)

Description	Trk B (Phospho-Tyr706/Tyr707) Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	Q16620
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of paraffin...