

Product Datasheet

ERK1/2 (Phospho-Thr202/Tyr204) Antibody (orb191866)

Description	ERK1/2 (Phospho-Thr202/Tyr204) Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	P27361
Expiration Date	12 months from date of receipt.