



Product Datasheet

APOH Antibody (orb192150)


Description	APOH Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	P02749
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formalin...



IHC-P of human stomach cancer tissue usi...



Western blot analysis of NB4 cell lysate...