

Product Datasheet

West Nile Virus (WNV) NS1 protein, His Tag (orb1921714)

Description	Flaviviruses belong to the family of flaviviridae, most of which are arboviruses mainly including Zi...
Tested Applications	CLIA, IA, LF
Form/Appearance	Purified,Liquid
Storage	For long term storage, the product should be stored at -20°C or lower.
Note	For research use only
Application notes	Recommended Applications: LFIA, CLIA, MIA
Sample Types	Test sample type: Whole Blood. Test sample source: Human
Purity	> 90% pure (SDS-PAGE).
Source	HEK293
MW	39.61 kDa
Expiration Date	6 months from date of receipt.