

## Product Datasheet

### Marburg GP1 (orb1921718)

<b>Description</b>	The Ebola virus is an enveloped, non-segmented, negative-strand RNA virus, belonging to the filoviri...
<b>Tested Applications</b>	CLIA, IA, LF
<b>Form/Appearance</b>	Purified, Liquid
<b>Storage</b>	For long term storage, the product should be stored at -20°C or lower.
<b>Note</b>	For research use only
<b>Sample Types</b>	Test sample type: Whole Blood. Test sample source: Human
<b>Purity</b>	> 90% pure (SDS-PAGE).
<b>Source</b>	HEK293
<b>MW</b>	26.02 kDa
<b>Expiration Date</b>	6 months from date of receipt.