

Product Datasheet

Ebola BDBV GP1/GP2 (orb1921719)

Description	The Ebola virus is an enveloped, non-segmented, negative-strand RNA virus, belonging to the filoviri...
Tested Applications	CLIA, IA, LF
Form/Appearance	Purified, Liquid
Storage	For long term storage, the product should be stored at -20°C or lower.
Note	For research use only
Sample Types	Test sample type: Whole Blood. Test sample source: Human
Purity	> 90% pure (SDS-PAGE).
Source	HEK293
MW	50.26 kDa
Expiration Date	6 months from date of receipt.