

Product Datasheet

BLLF1 g350 N-Terminal (orb1921729)

Description	As the most abundant envelop glycoprotein of EBV, gp350 is encoded by the BLLF1 gene, heavily glycos...
Tested Applications	CLIA, IA, LF
Form/Appearance	Purified, Liquid
Storage	For long term storage, the product should be stored at -20°C or lower.
Note	For research use only
Sample Types	Test sample type: Whole Blood. Test sample source: Human
Purity	> 90% pure (SDS-PAGE).
Source	HEK293
MW	52.25 kDa
Expiration Date	6 months from date of receipt.