

Product Datasheet

BLLF1 g350 ful length (orb1921730)

Description	As the most abundant envelop glycoprotein of EBV, gp350 is encoded by the BLLF1 gene, heavily glycos...
Tested Applications	CLIA, IA, LF
Form/Appearance	Purified, Liquid
Storage	For long term storage, the product should be stored at -20°C or lower.
Note	For research use only
Sample Types	Test sample type: Whole Blood. Test sample source: Human
Purity	> 90% pure (SDS-PAGE).
Source	HEK293
MW	89.13 kDa
Expiration Date	6 months from date of receipt.