

## Product Datasheet

# Herpes Simplex Virus (HSV)-2 Envelope Glycoprotein D (gD), His Tag (orb1921732)

<b>Description</b>	A total of 11 HSV envelope glycoproteins have been identified, some of them are strong inducers of a...
<b>Tested Applications</b>	CLIA, IA, LF
<b>Form/Appearance</b>	Purified,Liquid
<b>Storage</b>	For long term storage, the product should be stored at -20°C or lower.
<b>Note</b>	For research use only
<b>Application notes</b>	Recommended Applications: LFIA, CLIA, MIA
<b>Sample Types</b>	Test sample type: Whole Blood. Test sample source: Human
<b>Purity</b>	> 90% pure (SDS-PAGE).
<b>Source</b>	HEK293
<b>MW</b>	35.88kDa
<b>Expiration Date</b>	6 months from date of receipt.