

Product Datasheet

Herpes Simplex Virus (HSV)-1 Envelope Glycoprotein G1 (gG1), His Tag (orb1921733)

Description	A total of 11 HSV envelope glycoproteins have been identified, some of them are strong inducers of a...
Tested Applications	CLIA, IA, LF
Form/Appearance	Purified,Liquid
Storage	For long term storage, the product should be stored at -20°C or lower.
Note	For research use only
Application notes	Recommended Applications: LFIA, CLIA, MIA
Sample Types	Test sample type: Whole Blood. Test sample source: Human
Purity	> 90% pure (SDS-PAGE).
Source	HEK293
MW	17.39 kDa
Expiration Date	6 months from date of receipt.