

Product Datasheet

TpN15,17,47 (GST Tag) Protein (orb1921744)

Description	Treponema pallidum, a spirochete of the family Spirochaetaceae, is the etiologic agent of venereal s...
Tested Applications	CLIA, IA, LF
Form/Appearance	Purified,Liquid
Storage	For long term storage, the product should be stored at -20°C or lower.
Note	For research use only
Application notes	Recommended Applications: LFIA, CLIA, MIA. Test Methods: Self-produced antigen test; Test 10 Positive Reference Materials; Test Sample Extraction Solution, Test 10 normal Nasal Swab Samples
Sample Types	Test sample type: Serum / Plasma / Whole Blood. Test sample source: Human
Purity	> 90% pure (SDS-PAGE).
Source	E. coli
MW	35.83 kDa
Expiration Date	6 months from date of receipt.