

## Product Datasheet

### Adenovirus Monoclonal Antibody (orb1921796)

<b>Description</b>	Adenoviruses (AdVs) are non-enveloped icosahedral viruses with a double-stranded DNA genome. Hexon i...
<b>Species/Host</b>	Mouse
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	CLIA, LF
<b>Form/Appearance</b>	Purified,Liquid
<b>Storage</b>	For long term storage, the product should be stored at -20°C or lower.
<b>Note</b>	For research use only
<b>Application notes</b>	Recommended Applications: LFIA , CLIA. Test Methods: Self-produced antigen test; Test 10 Positive Reference Materials; Test Sample Extraction Solution, Test 10 normal Nasal Swab Samples; Test Covid-19, Influenza A/B & RSV controls
<b>Sample Types</b>	Test sample type: Nasopharyngeal swab /Oropharyngeal swab / Nasal swab / Serum / Plasma / Whole Blood / Stool. Test sample source: Human
<b>Clonality</b>	Monoclonal
<b>Purity</b>	> 95% pure (SDS-PAGE).
<b>MW</b>	1000 kDa
<b>Expiration Date</b>	12 months from date of receipt.