

Product Datasheet

UCP1 antibody (orb19218)

Description	Goat polyclonal antibody to UCP1
Species/Host	Goat
Reactivity	Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	UCP1
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:128000. WB: Approx 35kDa band observed in Human Adipose lysates (calculated MW of 33.0kDa according to NP_068605.1). Recommended concentration: 0.1-1 µg/ml.
Protein Sequence	EQLKRELSKSRQ
Clonality	Polyclonal
MW	33
Entrez	7350
Dilution Range	ELISA: 1:128000, WB: 0.5-1 µg/ml
Expiration Date	12 months from date of receipt.



orb19218
(1µg/ml)
staining of
Human Adip...