

## Product Datasheet

# Recombinant Human Respiratory Syncytial Virus (RSV) Recombinant F Protein (orb1921802)

<b>Description</b>	Respiratory syncytial virus (RSV) is a filamentous enveloped, negative-sense, single-stranded RNA vi...
<b>Tested Applications</b>	CLIA, LF
<b>Form/Appearance</b>	Purified,Liquid
<b>Storage</b>	For long term storage, the product should be stored at -20°C or lower.
<b>Note</b>	For research use only
<b>Application notes</b>	Recommended Applications: LFIA , CLIA. Test Methods: Test Covid-19, Influenza A/B & Aendovirus controls
<b>Sample Types</b>	Test sample type: Nasopharyngeal swab /Oropharyngeal swab / Nasal swab. Test sample source: Human
<b>Purity</b>	> 95% pure (SDS-PAGE).
<b>Source</b>	HEK293
<b>MW</b>	55.63 kDa
<b>Expiration Date</b>	6 months from date of receipt.