

Product Datasheet

Human Parainfluenza Virus 3 (HPIV3) Monoclonal Antibody (orb1921803)

Description	The human parainfluenza virus (HPIV) is an enveloped, negative-sense, single-stranded RNA virus that...
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	CLIA, LF
Form/Appearance	Purified,Liquid
Storage	For long term storage, the product should be stored at -20°C or lower.
Note	For research use only
Application notes	Recommended Applications: LFIA , CLIA. Test Methods: Self-produced antigen test; Test 10 Positive Reference Materials; Test Sample Extraction Solution, Test 10 normal Nasal Swab Samples; Test Covid-19, Influenza A/B, Aendovirus & RSV controls
Sample Types	Test sample type: Nasopharyngeal swab /Oropharyngeal swab / Nasal swab. Test sample source: Human
Clonality	Monoclonal
Purity	> 95% pure (SDS-PAGE).
MW	150 kDa
Expiration Date	12 months from date of receipt.