

## Product Datasheet

### Recombinant Human OSM(his Tag) (orb1921820)

<b>Description</b>	Oncostatin M (OSM) is a multifunctional cytokine that belongs to the Interleukin-6 subfamily. Among ...
<b>Endotoxins</b>	≤10 EU/mg
<b>Form/Appearance</b>	Lyophilized Powder
<b>Storage</b>	Use a manual defrost freezer and avoid repeated freeze - thaw cycles. 12 months from date of receipt, -20 to -70°C as supplied. 1 month, 2 to 8°C under sterile conditions after reconstitution. 6 months, -20 to -70°C under sterile conditions after reconstitution.
<b>Note</b>	For research use only
<b>Application notes</b>	Application Notes: For research use only. Minimum Sample Size: 50µg/vial
<b>Purity</b>	≥95%
<b>Source</b>	CHO Cell
<b>MW</b>	26.8 kDa
<b>Uniprot ID</b>	<a href="#">P13725</a>
<b>Entrez</b>	<a href="#">5008</a>
<b>Expiration Date</b>	6 months from date of receipt.