

Product Datasheet

Recombinant Human Vitronectin (orb1921827)

Description	Vitronectin (VN) is an abundant adhesive glycoprotein in blood plasma and is found associated with d...
Endotoxins	≤10 EU/mg
Form/Appearance	Lyophilized Powder
Storage	Use a manual defrost freezer and avoid repeated freeze - thaw cycles. 12 months from date of receipt, -20 to -70°C as supplied. 1 month, 2 to 8°C under sterile conditions after reconstitution. 6 months, -20 to -100°C under sterile conditions after reconstitution.
Note	For research use only
Application notes	Application Notes: For research use only. Minimum Sample Size: 1 mg/vial
Purity	≥90%
Source	CHO Cell
MW	52.10 kDa
Uniprot ID	P04004
Entrez	7448
Expiration Date	6 months from date of receipt.