

Product Datasheet

Recombinant Human TGF- β 1 (orb1921828)

Description	Transforming growth factor- β 1 (TGF- β 1) was the first fully cloned of the TGF- β superfamily member...
Endotoxins	\leq 10 EU/mg
Form/Appearance	Lyophilized Powder
Storage	Use a manual defrost freezer and avoid repeated freeze - thaw cycles. 12 months from date of receipt, -20 to -70°C as supplied. 1 month, 2 to 8°C under sterile conditions after reconstitution. 6 months, -20 to -97°C under sterile conditions after reconstitution.
Note	For research use only
Application notes	Application Notes: For research use only. Minimum Sample Size: 50 μ g/vial
Purity	\geq 95%
Source	CHO Cell
MW	12.62 kDa
Uniprot ID	P01137
Entrez	7040
Expiration Date	6 months from date of receipt.