

Product Datasheet

Recombinant Human LR3 IGF-1 (orb1921829)

Description	Long R3 IGF-I (LR3 IGF-I) is a 9.2 kDa synthetic analog of IGF-I that is generated by modifying the ...
Endotoxins	≤10 EU/mg
Form/Appearance	Lyophilized Powder
Storage	Use a manual defrost freezer and avoid repeated freeze - thaw cycles. 12 months from date of receipt, -20 to -70°C as supplied. 1 month, 2 to 8°C under sterile conditions after reconstitution. 6 months, -20 to -94°C under sterile conditions after reconstitution.
Note	For research use only
Application notes	Application Notes: For research use only. Minimum Sample Size: 1 mg/vial
Purity	≥95%
Source	CHO Cell
MW	9kDa
Uniprot ID	P05019
Entrez	3479
Expiration Date	6 months from date of receipt.