

Product Datasheet

Recombinant Human VEGF165 (His Tag) (orb1921837)

Description	VEGF165 Protein, Human (P.pastoris) is a Vascular Endothelial Growth Factor A (VEGF-A) isoform. VEGF...
Endotoxins	≤10 EU/mg
Form/Appearance	Lyophilized Powder
Storage	Use a manual defrost freezer and avoid repeated freeze - thaw cycles. 12 months from date of receipt, -20 to -70°C as supplied. 1 month, 2 to 8°C under sterile conditions after reconstitution. 6 months, -20 to -86°C under sterile conditions after reconstitution.
Note	For research use only
Application notes	Application Notes: For research use only. Minimum Sample Size: 50µg/vial
Purity	≥95%
Source	CHO Cell
MW	20.1kDa
Uniprot ID	P15692
Entrez	7422
Expiration Date	6 months from date of receipt.