

Product Datasheet

Recombinant Human Fibronectin (orb1921841)

Description	Fibronectin, also known as fibronectin, English name fibronectin (abbreviated FN), is a polymer glyc...
Endotoxins	≤100 EU/mg
Form/Appearance	Lyophilized Powder
Storage	Use a manual defrost freezer and avoid repeated freeze - thaw cycles. 12 months from date of receipt, -20 to -70°C as supplied. 1 month, 2 to 8°C under sterile conditions after reconstitution. 6 months, -20 to -81°C under sterile conditions after reconstitution.
Note	For research use only
Application notes	Application Notes: For research use only. Minimum Sample Size: 1 mg/vial
Purity	≥95%
Source	E.coli
MW	62.6kDa
Uniprot ID	P02751
Entrez	2335
Expiration Date	6 months from date of receipt.