

Product Datasheet

Recombinant Human IL-21(hFc Tag) (orb1921843)

Description	Interleukin (IL) -21 is a member of the common gamma chain (γ c) family of cytokines and is expresse...
Endotoxins	≤ 10 EU/mg
Form/Appearance	Lyophilized Powder
Storage	Use a manual defrost freezer and avoid repeated freeze - thaw cycles. 12 months from date of receipt, -20 to -70°C as supplied. 1 month, 2 to 8°C under sterile conditions after reconstitution. 6 months, -20 to -79°C under sterile conditions after reconstitution.
Note	For research use only
Application notes	Application Notes: For research use only. Minimum Sample Size: 50 μ g/vial
Purity	$\geq 90\%$
Source	CHO Cell
MW	83.77 kDa
Uniprot ID	Q9HBE4
Entrez	59067
Expiration Date	6 months from date of receipt.