

## Product Datasheet

### Recombinant Human IL-21 (orb1921844)

<b>Description</b>	Interleukin (IL) -21 is a member of the common gamma chain ( $\gamma$ c) family of cytokines and is expresse...
<b>Endotoxins</b>	$\leq 10$ EU/mg
<b>Form/Appearance</b>	Lyophilized Powder
<b>Storage</b>	Use a manual defrost freezer and avoid repeated freeze - thaw cycles. 12 months from date of receipt, $-20$ to $-70^{\circ}\text{C}$ as supplied. 1 month, $2$ to $8^{\circ}\text{C}$ under sterile conditions after reconstitution. 6 months, $-20$ to $-79^{\circ}\text{C}$ under sterile conditions after reconstitution.
<b>Note</b>	For research use only
<b>Application notes</b>	Application Notes: For research use only. Minimum Sample Size: $50\mu\text{g/vial}$
<b>Purity</b>	$\geq 95\%$
<b>Source</b>	CHO Cell
<b>MW</b>	83.77 kDa
<b>Uniprot ID</b>	<a href="#">Q9HBE4</a>
<b>Entrez</b>	<a href="#">59067</a>
<b>Expiration Date</b>	6 months from date of receipt.