

## Product Datasheet

### Recombinant Human IL-15(hFc Tag) (orb1921846)

<b>Description</b>	Interleukin (IL)-15 is a 14-15 kDa four-helix bundle cytokine. IL-15, unlike other cytokines, is wid...
<b>Endotoxins</b>	≤10 EU/mg
<b>Form/Appearance</b>	Lyophilized Powder
<b>Storage</b>	Use a manual defrost freezer and avoid repeated freeze - thaw cycles. 12 months from date of receipt, -20 to -70°C as supplied. 1 month, 2 to 8°C under sterile conditions after reconstitution. 6 months, -20 to -77°C under sterile conditions after reconstitution.
<b>Note</b>	For research use only
<b>Application notes</b>	Application Notes: For research use only. Minimum Sample Size: 50µg/vial
<b>Purity</b>	≥90%
<b>Source</b>	CHO Cell
<b>MW</b>	77.33 kDa
<b>Uniprot ID</b>	<b>P40933</b>
<b>Entrez</b>	<b>3600</b>
<b>Expiration Date</b>	6 months from date of receipt.