

Product Datasheet

TAGLN/Transgelin Rabbit mAb Antibody (orb1921965)

Description	TAGLN/Transgelin Rabbit mAb Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	A synthesized peptide derived from human TAGLN
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.163.
Form/Appearance	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Uniprot ID	Q01995
Dilution Range	IHC: 1/100-1/200
Expiration Date	12 months from date of receipt.