

## Product Datasheet

### LI Cadherin Rabbit mAb Antibody (orb1922315)

<b>Description</b>	LI Cadherin Rabbit mAb Antibody
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Human
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	IHC-P
<b>Immunogen</b>	A synthesized peptide derived from human LI Cadherin.
<b>Preservatives</b>	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.
<b>Form/Appearance</b>	Liquid
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Isotype</b>	IgG
<b>Clonality</b>	Monoclonal
<b>Uniprot ID</b>	<a href="#">Q12864</a>
<b>Dilution Range</b>	IHC: 1/100-1/200
<b>Expiration Date</b>	12 months from date of receipt.