



Product Datasheet

CD47 Antibody (orb192374)


Description	CD47 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Concentration	1.0mg/ml
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	Q08722
Dilution Range	WB 1:500 - 1:2000, IHC 1:50 - 1:200, IF 1:50 - 1:200
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of U-2 OS ce...



Immunofluorescence analysis of L929 cell...



Western blot analysis of extracts of var...