

## Product Datasheet

### Transgelin Rabbit mAb Antibody (orb1923834)

|                            |  |
|----------------------------|--|
| <b>Description</b>         | Transgelin Rabbit mAb Antibody   |
| <b>Species/Host</b>        | Rabbit   |
| <b>Reactivity</b>          | Human  |
| <b>Conjugation</b>         | Unconjugated   |
| <b>Tested Applications</b> | IHC-P, WB  |
| <b>Immunogen</b>           | A synthesized peptide derived from human TAGLN   |
| <b>Preservatives</b>       | Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. |
| <b>Form/Appearance</b>     | Liquid   |
| <b>Storage</b>             | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.           |
| <b>Note</b>                | For research use only  |
| <b>Isotype</b>             | IgG  |
| <b>Clonality</b>           | Monoclonal   |
| <b>MW</b>                  | Calculated MW: 23 kDa; Observed MW: 23 kDa   |
| <b>Uniprot ID</b>          | <a href="#">Q01995</a>   |
| <b>Dilution Range</b>      | WB: 1/500-1/1000 IHC: 1/50-1/100   |
| <b>Expiration Date</b>     | 12 months from date of receipt.  |