

Product Datasheet

Retinoid X Receptor alpha Rabbit mAb Antibody (orb1923898)

Description	Retinoid X Receptor alpha Rabbit mAb Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IP, WB
Immunogen	A synthesized peptide derived from human Retinoid X Receptor alpha
Preservatives	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
MW	Calculated MW: 51 kDa; Observed MW: 51 kDa
Uniprot ID	P19793
Dilution Range	WB: 1/500-1/1000 IF: 1/50-1/200 IP: 1/50
Expiration Date	12 months from date of receipt.