

## Product Datasheet

### Goat Anti-Rabbit IgG H&L (Cy3) (orb1923988)

<b>Description</b>	Goat Anti-Rabbit IgG H&L (Cy3)
<b>Species/Host</b>	Goat
<b>Reactivity</b>	Rabbit
<b>Conjugation</b>	Cy3
<b>Tested Applications</b>	ELISA, FC, ICC, IHC-Fr, IHC-P
<b>Immunogen</b>	Rabbit IgG
<b>Preservatives</b>	Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% rAlbumin, 50% glycerol, 0.02% Sodium Azide
<b>Form/Appearance</b>	Liquid
<b>Concentration</b>	1 mg/ml
<b>Storage</b>	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
<b>Note</b>	For research use only
<b>Clonality</b>	Polyclonal
<b>Dilution Range</b>	IF: 1/100-1/1000 FC: 1/500-1/1000 ELISA: 1/5000-1/20000
<b>Expiration Date</b>	12 months from date of receipt.