

Product Datasheet

Goat Anti-Mouse IgG H&L (Cy3) (orb1923989)

Description	Goat Anti-Mouse IgG H&L (Cy3)
Species/Host	Goat
Reactivity	Mouse
Conjugation	Cy3
Tested Applications	ELISA, FC, ICC, IHC-Fr
Immunogen	Mouse IgG
Preservatives	Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% rAlbumin, 50% glycerol, 0.02% Sodium Azide
Form/Appearance	Liquid
Concentration	1 mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Dilution Range	IF: 1/100-1/1000 FC: 1/500-1/1000 ELISA: 1/5000-1/20000
Expiration Date	12 months from date of receipt.