

Product Datasheet

Rabbit Anti-Rat IgG H&L (FITC) (orb1923990)

Description	Rabbit Anti-Rat IgG H&L (FITC)
Species/Host	Rabbit
Reactivity	Rat
Conjugation	FITC
Tested Applications	FC, ICC
Immunogen	Rat IgG
Preservatives	Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% rAlbumin, 50% glycerol, 0.05% Sodium Azide
Form/Appearance	Liquid
Concentration	1 mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Dilution Range	IF: 1/100-1/1000 FC: 1/500-1/1000
Expiration Date	12 months from date of receipt.