

Product Datasheet

Goat Anti-Rabbit IgG H&L (Biotin) (orb1923997)

Description	Goat Anti-Rabbit IgG H&L (Biotin)
Species/Host	Goat
Reactivity	Rabbit
Conjugation	Biotin
Tested Applications	ELISA, ICC, IHC-Fr, IHC-P, WB
Immunogen	Rabbit IgG
Preservatives	Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% rAlbumin, 50% glycerol, 0.05% Sodium Azide
Form/Appearance	Liquid
Concentration	1 mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB: 1/5000-1/20000 IF: 1/100-1/1000 ELISA: 1/5000-1/20000
Expiration Date	12 months from date of receipt.