

Product Datasheet

Hepatitis B Surface Antigen Rabbit mAb (orb1924037)

Description	Hepatitis B Surface Antigen Rabbit mAb
Species/Host	Rabbit
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	A synthesized peptide derived from Hepatitis B Surface Antigen
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Uniprot ID	P03141
Dilution Range	IHC: 1/100-1/200
Expiration Date	12 months from date of receipt.