



Product Datasheet

Leptin Receptor antibody (orb19243)

Description	Goat polyclonal antibody to Leptin Receptor
Species/Host	Goat
Reactivity	Canine, Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IHC, WB
Target	Leptin Receptor
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:16000. WB: Approx 150kDa band observed in Human Brain lysates (calculated MW of 132kDa according to NP_002294.2). Recommended concentration: 1-3 µg/ml.
Protein Sequence	TQDDIEKHQSDAG
Clonality	Polyclonal
MW	132; 102; 109
Entrez	3953
Dilution Range	ELISA: 1:16000, WB: 0.1-1 µg/ml, IHC-P: 4 µg/ml
Expiration Date	12 months from date of receipt.



orb19243
Flow
cytometric
analysis of
par...



orb19243
Negative
Control
showing
staini...



orb19243
(6µg/ml)
staining of
paraffin ...