

Product Datasheet

[KO] GIPC1 Mouse mAb Antibody (orb1924327)

Description	[KO] GIPC1 Mouse mAb Antibody
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Full length human recombinant protein of human GIPC1 produced in E.coli.
Preservatives	Liquid in PBS containing 50% glycerol, 1% rAlbumin and 0.02% sodium azide, pH 7.3.
Form/Appearance	Liquid
Concentration	1 mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
MW	Calculated MW: 36 kDa; Observed MW: 42 kDa
Uniprot ID	O14908
Dilution Range	WB: 1/500-1/1000 IHC: 1/50-1/100
Expiration Date	12 months from date of receipt.