


## Product Datasheet

### PRNP antibody (orb19246)

<b>Description</b>	Goat polyclonal antibody to PRNP
<b>Species/Host</b>	Goat
<b>Reactivity</b>	Bovine, Human, Porcine
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA, WB
<b>Target</b>	Prion Protein (143-153)
<b>Preservatives</b>	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
<b>Storage</b>	Aliquot and store at -20°C. Minimize freezing and thawing.
<b>Note</b>	For research use only
<b>Application notes</b>	ELISA: Peptide ELISA: antibody detection limit dilution 1:64000.WB: Approx 25kDa band observed in Human Brain lysates (calculated MW of 27.7kDa according to NP_000302.1 and NP_898902.1). Recommended concentration: 0.3-1 µg/ml.
<b>Protein Sequence</b>	SDYEDRYREN
<b>Clonality</b>	Polyclonal
<b>MW</b>	27.7
<b>Entrez</b>	<a href="#">5621</a>
<b>Dilution Range</b>	ELISA: 1:64000, WB: 0.3-1 µg/ml
<b>Expiration Date</b>	12 months from date of receipt.



orb19246  
(0.3ug/ml)  
staining of  
Human Br...