

Product Datasheet

Goat Anti-Rabbit IgG H&L \square HRP \square (orb1924910)

Description	Goat Anti-Rabbit IgG H&L \square HRP \square
Species/Host	Goat
Reactivity	Rabbit
Conjugation	HRP
Tested Applications	ELISA, IHC-P, WB
Immunogen	Rabbit IgG
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.
Form/Appearance	liquid
Concentration	1 mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB: 1/5000-1/20000 IHC: 1/500-1/1000 ELISA: 1/5000-1/20000
Expiration Date	12 months from date of receipt.