

Product Datasheet

Goat Anti-Rabbit IgG H&L (FITC) (orb1924912)

Description	Goat Anti-Rabbit IgG H&L (FITC)
Species/Host	Goat
Reactivity	Rabbit
Conjugation	FITC
Tested Applications	FC, ICC
Immunogen	Rabbit IgG
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.
Form/Appearance	Liquid
Concentration	1 mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Dilution Range	IF: 1/100-1/1000 FC: 1/500-1/1000
Expiration Date	12 months from date of receipt.