

Product Datasheet

HA Tag Mouse mAb (HRP) Antibody (orb1924920)

Description	HA Tag Mouse mAb (HRP) Antibody
Species/Host	Mouse
Reactivity	Other
Conjugation	HRP
Tested Applications	ELISA, WB
Immunogen	YPYDVPDYA conjugated to KLH.
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.
Form/Appearance	Liquid
Concentration	1 mg/ml
Storage	Store at 4°C for at least 6 months. Avoid freeze/thaw cycles
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Dilution Range	WB: 1/1000-1/10000 ELISA: 1/5000-1/20000
Expiration Date	12 months from date of receipt.