

Product Datasheet

HA Tag (9A3) Mouse mAb Antibody (orb1924922)

Description	HA Tag (9A3) Mouse mAb Antibody
Species/Host	Mouse
Reactivity	Other
Conjugation	Unconjugated
Tested Applications	ICC, IP, WB
Immunogen	Synthetic Peptide of HA-Tag
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.
Form/Appearance	Liquid
Concentration	1 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Dilution Range	WB: 1/1000-1/10000 IF: 1/100-1/500 IP: 1/200
Expiration Date	12 months from date of receipt.