

Product Datasheet

c-Myc Tag (9C4) Mouse mAb Antibody (orb1924933)

Description	c-Myc Tag (9C4) Mouse mAb Antibody
Species/Host	Mouse
Reactivity	Other
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	Monoclonal antibodies are produced by immunizing mouse with synthetic peptide (EQKLISEEDL) coupled to KLH.
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.
Form/Appearance	Liquid
Concentration	0.5 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Dilution Range	WB: 1/1000-1/10000 IP: 1/500
Expiration Date	12 months from date of receipt.