

Product Datasheet

c-Myc Tag (9E10) Mouse mAb (HRP) Antibody (orb1925038)

Description	c-Myc Tag (9E10) Mouse mAb (HRP) Antibody
Species/Host	Mouse
Reactivity	Other
Conjugation	HRP
Tested Applications	ELISA, IHC-P, IP, WB
Immunogen	Synthetic peptide conjugated to Keyhole Limpet Haemocyanin (KLH).
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.
Form/Appearance	Liquid
Concentration	1 mg/ml
Storage	Store at 4°C for at least 6 months. Avoid freeze/thaw cycles
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Dilution Range	WB: 1/1000-1/10000 IHC: 1/50-1/200 IP: 1/500 ELISA: 1/5000-1/20000
Expiration Date	12 months from date of receipt.