

Product Datasheet

Deamidated gliadin peptide antigen (orb1925045)

Description	This recombinant protein contains multiple deamidated peptide of gliadin with a calculated MW of 37....
Tested Applications	CLIA, ELISA, LF, WB
Preservatives	0.1% ProClin300. Tris-buffered saline, pH 8.0
Form/Appearance	Purified, Liquid
Storage	Upon receipt, quick spin before open vial, aliquot and store at $\leq -20^{\circ}\text{C}$. Avoid multiple freeze/thaw cycles
Note	For research use only
Purity	> 90% pure by SDS-PAGE
Expiration Date	6 months from date of receipt.